

Heelift® – heel off-loading 3 versions for the best treatment



Effective off-loading for heel pressure ulcers, especially for bedridden patients

Heelift®



> The Original model – an ideal and economical device for heel off-loading

Heelift® Glide



> The premium version of the Heelift – providing an additional forefoot strap for more support and a gliding exterior

Heelift® AFO



> Heelift AFO – improved version combining all features of the Heelift Glide with a stabilizing plate for limited weight bearing

2 types of interior foam available



Standard
convoluted foam for
better air circulation




Smooth
smooth foam for
sensitive skin

- > for the prevention and treatment of heel pressure ulcers, diabetic foot syndrome and decubitus
- > diabetes with neurological sensory deficit
- > prolonged immobility or coma
- > allows for patient-specific customization



DARCO

What makes the Heelift® the right boot for your patient?



All Heelift® products

- 1 Ventilation Holes provide additional air circulation to promote increased patient comfort.
- 2 Heel Suspension Opening allows for visual inspection of the heel.
- 3 Affixed Elevation Pad lifts the leg to suspend the heel, promoting the ideal environment for an ulcer to heal.
- 4 Soft »easy-on, easy-off« Closure Straps are a simple solution for a custom fit.
- 5 Friction-free Tricot Backing aids in patient mobility and helps ensure patient safety.

Heelift® Glide

- 6 Durable Exterior allows the patient's leg to glide easily across bed sheets.
- 7 Forefoot Strap keeps the boot in place and helps support the forefoot.

Heelift® AFO

- 6 Durable Exterior allows the patient's leg to glide easily across bed sheets.
- 7 Forefoot Strap keeps the boot in place and helps support the forefoot.
- 8 Semi-rigid brace provides the best possible control of the ankle's position and motion.



Standard



Smooth

Choice of two interiors

The **standard** version with convoluted foam material allows better air circulation. The attached nubs provide a pleasant cooling effect.

The **smooth** foam material was developed especially for patients with sensitive skin. This helps prevent the formation of oedema. This version is recommended for diabetic patients.

SIZE	CIRCUMFERENCE CALF	ART. NO.	
Heelift®	Color: Beige	Standard	Smooth
Petite	15 – 25 cm	DM-HLP	DM-HLSP
Normal	20 – 36 cm	DM-HL	DM-HLS
Bariatric	30 – 58 cm	–	DM-HLSB
Heelift® Glide	Color: Red	Standard	Smooth
Petite	15 – 25 cm	DM-HLPG	DM-HLSPG
Normal	20 – 36 cm	DM-HLG	DM-HLSG
Bariatric	30 – 58 cm	–	DM-HLSBG
Heelift® AFO	Color: Blue	Standard	Smooth
Universal	20 – 36 cm	DM-HLAFO	DM-HLSAFO

All Heelift products fit right and left. They can be washed safely at 50° C and can be autoclaved for sterility.