

FAST DRY

derma+flex[®]

TOPICAL SKIN ADHESIVE

A High Viscosity Tissue Adhesive Like No Other!

FAST DRYING 60/40 BLENDED FORMULATION

- Optimal Flexibility
- No Cracking or Breaking

MICROBIAL BARRIER

- Protects the affected area by preventing bacteria from penetrating the protective barrier, greatly reducing the likelihood of infection*

QUICK DRYING

- Polymerizes (dries) quickly, usually within 60 seconds

PAINLESS APPLICATION:

- No sting or burn

WATERPROOF FORMULA

- Prevents water from reaching the wound
- Allows your patient to shower within hours of their surgery

HIGH VISCOSITY

- Syrup-like Viscosity at ~ 250cps
- Does not run or drip when applied correctly
- The perfect choice for vertical surface applications

MULTIPLE USES ON A SINGLE PATIENT

- Full size product (0.7 g)

IMPROVED DUAL FUNCTION APPLICATOR

- Tilt on edge for precise lines
- Tilt down for broad coverage

* As long as adhesive is kept intact



A product of Chemence Medical

www.dermaflex.com